SHOULD WE MEASURE OR CALCULATE FREE 250H VITAMIN D IN VITAMIN D RESEARCH?

N. Heureux^{*} DIAsource Immunoassays, Louvain-la-Neuve, Belgium

Free 25OH Vitamin D represents 0.04% of the total 25OH Vitamin D


Almost all circulating 250H Vitamin D in serum is bound to Vitamin D Binding Protein DBP (88%) and Albumin (12%). A very small fraction, approximately 0.04 % of the 25OH Vitamin D, circulates in the free, non-protein bound form.



× Context: compare direct quantitation to calculated levels in cirrhotic patients, pregnant women (2nd and 3^d trimester) and control adults.


× Methods: total 250HD by LC-MS, Albumin by BCG and DBP by ELISA (MoAb).

× Results: calculations overestimated free 250H D levels compared to directly measured with larger fold differences in African Americans compared to Caucasians. Directly measured free 250HD was correlated with iPTH and Ca, but calculated was not.



	ELISA	Calc.
iPTH	p<0.02	p = 0.46
Ca	p<0.004	p = 0.47

Schwartz J.B. et al. J. Clin. Endocrinol. Metab. 2014, 99(5), 1631-7.



Conditions affecting the free 250H Vitamin D concentration

The concentration of DBP is not constant and can be influenced by a number of factors including

Obesity

Pregnancy (\times)

× The use of oral contraceptives

Liver disease

Renal disease

Proteinuria

2016 – Comparison in African **Americans and Caucasians**

Context: free 250HD was measured by ELISA and calculated (X) amongst ethnical groups.

× Methods: total 250HD by LC-MS, Albumin by BCG and DBP by polyclonal radial immunodiffusion (pRID) and ELISA (MoAb and PoAb).

Results: calculated free 25OHD concentrations were higher than directly measured values. Since calculated free 250HD derives from measures of DBP, the accuracy of DBP assays is critical but not constant.

× Hormone replacement therapy **(** ×)



In case of elevated concentration of DBP the % of free 25OH Vitamin D is decreased. In case of decreased concentration of DBP the % of free 25OH Vitamin D is increased.

(×) The supplementation with Vitamin D2 or D3 also affects differently the free 250H Vitamin D levels.


Methods for the evaluation of free 250H Vitamin D

⊗ Centrifugal ultrafiltration

- Time and labor consuming
- Uses H³ and C¹⁴
- Uses a few assumptions

Calculations

- Requires 3 assays per sample (DBP, Albumin, total 250HD)



Orwoll E.S. et al. J. Clin. Endocrinol. Metab. 2016, 101(5), 2226-34.



Context: seasonal variation of 250HD was evaluated by (X) direct measurement and by calculations.

(×) Methods: total 250HD by CLIA, Albumin by Cobas and DBP by ELISA (MoAb). × Results: the directly measured free 250HD was significantly lower than the calculated, and directly measured showed a stronger correlation with serum iPTH (P<0.001) than calculated.

future diagnostics

- Should ideally include DBP genotyping
- DBP assays are not of equal quality

✓ Direct assay - ELISA

- Direct measurement method
- Accessible to any laboratory
- Can be automated for large cohorts

Klingberg E. et al. Endocr. Connect. 2017, 6(2), 111-120.

Directly measured free 250HD is lower than calculated Directly measured free 250HD is better correlated to PTH and Ca

D. Source

Join the Scientific Community who talks about free Vitamin D: LinkedIn group "Free Vitamin D"