

FREE 25OH VITAMIN D ELISA

> SIMPLE-ACCURATE-AND-DIRECT

→ WHY YOU SHOULD USE THE DIASOURCE FREE VITAMIN D KIT

- ① According to the free hormone hypothesis the biological activity of the hormone is directly linked to the concentration of its free form
 - 25OH Vitamin D is for >99.9% bound to binding proteins: 90% to DBP, 10% to Albumin
 - About 0.04% circulates as the free form = "free 25OH Vitamin D"
 - All current 25OH Vitamin D assays measure the sum of the bound and free forms
- ② Free 25OH Vitamin D seems to be a better marker of Vitamin D status than total 25OH Vitamin D for:
 - Conditions affecting the binding proteins concentrations
 - Obesity/Insulin
 - Pregnancy
 - Cancer
 - Respiratory disease
 - Liver disease
 - Renal disease
 - Osteoporosis/Bone mineral density
 - Intensive care
 - Ethnic groups having polymorphic forms of DBP with different affinities for 25OH Vitamin D
 - Black
 - Hispanic
 - Asian
- ⊗ Following methods for measuring free 25OH Vitamin D exist:
 - Centrifugal ultrafiltration – accurate but long and tedious
 - Calculations – requires 3 assays and can be inaccurate
 - Direct measurement by ELISA – simple and accurate
- ⊗ DIAsource offers the only direct measurement assay in the market as an RUO, but
 - CE-MARKED
 - FDA REGISTRATION EXPECTED IN 2018

> OUR OFFERING

If you need support for internal validation and/or for a study design of our free Vitamin D ELISA test in your laboratory please contact Diasource:

+ COMMERCIAL INFORMATION

Commercial: Mr. Jan Wauters
+32 (0)479 94 34 78
jan.wauters@diasource.be

Business Development: Aziza El-Bouyahyaoui
+33 (0)6 37 19 84 52
aziza.elmahdaoui@diasource.be

+ SCIENTIFIC OR TECHNICAL QUESTION?

Scientific or technical: Dr. Nicolas Heureux
+32 (0)10 84 99 40
nicolas.heureux@diasource.be

+ ORDERING INFORMATION

Contact our Customer Service
Tel: +32 (0)10 84 99 00
Fax: +32 (0)10 84 99 96
Email: customer.service@diasource.be

CAT#	Product Description	Format
KAPF1991	Free 25OH Vitamin D ELISA	Kit 96 assays

⊗ REFERENCES:

- BOUILLON R. (2016), Free 25-hydroxyvitamin D: impact of vitamin D binding protein assays on racial genotypic associations, JCEM, <http://dx.doi.org/10.1210/jc.2016-1104>.
- SOLLID S.T. (2016), Effects of vitamin D binding protein phenotypes and vitamin D supplementation on serum total 25(OH) D and directly measured free 25(OH)D, Eur. J. Endocrinol. April 1, 174:445-452.
- TANGPRICHA V. (2015), Free 25-Hydroxyvitamin D Concentrations in Cystic Fibrosis, Am. J. Med. Sci. 2015 Nov;350(5):374-9.
- ALOIA J. (2015), Free 25(OH)D and the Vitamin D Paradox in African Americans, J. Clin. Endocrinol. Metab. 2015 Jul 10;JC20152066.
- SCHWARTZ J.B. (2014), A comparison of direct and calculated free 25(OH) Vitamin D levels in clinical populations, J. Clin. Endocrinol. Metab., 99(5):1631-7.
- BIKLE D. (2013), Variability in free 25(OH) vitamin D levels in clinical populations, J. Steroid Biochem. Mol. Biol., S0960-0760.

More references in our Technical Review 2014-02.

Manufactured by



Exclusively Distributed by

